

TECHNICAL DATA SHEET

CyFlow™ TCR gamma/delta PE Anti-Rt; Clone V65

REF AC827019

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Specifications

Antigen	TCR γ/δ
Alternative Names	—
Clone	V65
Clonality	monoclonal
Format	PE
Host / Isotype	Mouse / IgG1
Species Reactivity	Rat
Negative Species Reactivity	—
Quantity [Concentration]	0.1 mg [0.5 mg/ml]
Immunogen	TCR α/β -negative CD3-positive rat T cell hybridoma III.89.1.4 line

Specificity

The mouse monoclonal antibody V65 recognizes rat TCR γ/δ antigen, the subtype of T cell receptor expressed mainly in epithelial tissues and at the sites of infection.

Contact Information:

Sysmex Partec GmbH • Am Flugplatz 13 • 02828 Görlitz • Germany
Tel +49 3581 8746 0 • Fax +49 3581 8746 70 • E-mail: info@sysmex-partec.com

Application

The reagent is designed for Flow Cytometry analysis. Working concentrations should be determined by the investigator.

Other usages may be determined from the scientific literature.

Storage Buffer

The reagent is provided in phosphate buffered saline (PBS) solution, pH \approx 7.4, containing 0.1% (w/v) sodium azide.

Storage and Stability

Storage	Avoid prolonged exposure to light. Store in the dark at 2-8°C. Do not freeze.
Stability	Do not use after expiration date stamped on vial label.

Background Information

The antigen-specific T cell receptor (TCR) is composed of either α and β subunit, or γ and δ subunit. Majority of T cells present in the blood, lymph and secondary lymphoid organs express TCR α/β heterodimers, whereas the T cells expressing TCR γ/δ heterodimers are localized mainly in epithelial tissues and at the sites of infection. The subunits of TCR heterodimers are covalently bonded and in the endoplasmic reticulum they associate with CD3 subunits to form functional TCR-CD3 complex. Lack of expression of any of the chains is sufficient to stop cell surface expression.

References

- Giezeman-Smiths KM, Jonges LE, Chambers WH, Brisette-Storkus CS, Van Vlierberghe RL, Van Eendenburg JD, Eggermont AM, Fleuren GJ, Kuppen PJ: Novel monoclonal antibodies against membrane structures that are preferentially expressed on IL-2-activated rat NK cells. J Leukoc Biol. 1998 Feb; 63(2):209-15. < PMID: 9468279 >
- Tam FW, Smith J, Morel D, Karkar AM, Thompson EM, Cook HT, Pusey CD: Development of scarring and renal failure in a rat model of crescentic glomerulonephritis. Nephrol Dial Transplant. 1999 Jul; 14(7):1658-66. < PMID: 10435873 >
- Hosea HJ, Rector ES, Taylor CG: Zinc-deficient rats have fewer recent thymic emigrant (CD90+) T lymphocytes in spleen and blood. J Nutr. 2003 Dec; 133(12):4239-42. < PMID: 14652378 >

Contact Information:

Sysmex Partec GmbH • Am Flugplatz 13 • 02828 Görlitz • Germany
Tel +49 3581 8746 0 • Fax +49 3581 8746 70 • E-mail: info@sysmex-partec.com

- Henriksson G, Helgeland L, Midtvedt T, Stierna P, Brandtzaeg P: Immune response to *Mycoplasma pulmonis* in nasal mucosa is modulated by the normal microbiota. *Am J Respir Cell Mol Biol.* 2004 Dec; 31(6):657-62. < PMID: 15347559 >
- Xystrakis E, Dejean AS, Bernard I, Druet P, Liblau R, Gonzalez-Dunia D, Saoudi A: Identification of a novel natural regulatory CD8 T-cell subset and analysis of its mechanism of regulation. *Blood.* 2004 Nov 15; 104(10):3294-301. < PMID: 15271801 >
- van den Brandt J, Kwon SH, Hünig T, McPherson KG, Reichardt HM: Sustained pre-TCR expression in Notch1IC-transgenic rats impairs T cell maturation and selection. *J Immunol.* 2005 Jun 15; 174(12):7845-52. < PMID: 15944289 >
- Moghaddami M, Cleland LG, Mayrhofer G: MHC II+ CD45+ cells from synovium-rich tissues of normal rats: phenotype, comparison with macrophage and dendritic cell lineages and differentiation into mature dendritic cells in vitro. *Int Immunol.* 2005 Aug; 17(8):1103-15. < PMID: 16030130 >
- De Bernardis F, Lucciarini R, Boccanera M, Amantini C, Arancia S, Morrone S, Mosca M, Cassone A, Santoni G: Phenotypic and functional characterization of vaginal dendritic cells in a rat model of *Candida albicans* vaginitis. *Infect Immun.* 2006 Jul; 74(7):4282-94. < PMID: 16790803 >

The Safety Data Sheet for this product is available at www.sysmex-partec.com/services.

Contact Information:

Sysmex Partec GmbH • Am Flugplatz 13 • 02828 Görlitz • Germany
Tel +49 3581 8746 0 • Fax +49 3581 8746 70 • E-mail: info@sysmex-partec.com